

Product Datasheet

Anti-phospho-HSP70 (Tyr611)/PE-Cy5.5 (orb938426)

Description	Anti-phospho-HSP70 (Tyr611)/PE-Cy5.5
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	PE/Cy5.5
Tested Applications	FC, ICC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human HSP70 around the phosphorylation site of Tyr611
Target	phospho-HSP70 (Tyr611)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	PODMV8
Dilution Range	Flow-Cyt=1µg/Test,ICC=1:100
Expiration Date	12 months from date of receipt.