

Product Datasheet

NMDAR2B (phospho-Tyr1472) antibody (FITC) [Out of stock] (orb9387)

Description

Rabbit polyclonal to NMDAR2B

	conjugated to FITC
Species/Host	Rabbit
Conjugation	FITC
Tested Applications	FC, IF
Immunogen	KLH conjugated synthetic peptide derived from mouse NMDAR2B
Target	NMDAR2B (phospho-Tyr1472)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid: 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified by Protein A
MW	163 kDa
CAS Number	Glutamate receptor ionotropic, NMDA 2B
Uniprot ID	Q13224
Dilution Range	Flow-Cyt: 1:50-200, IF: 1:50-200
Expiration Date	12 months from date of receipt.