

Product Datasheet

DpnI (orb93993)

Description	Neoschizomers: BfuCI, Bsp143I, BstENII, BstKTI, BstMBI, DpnII, Kzo9I, NdeI
Form/Appearance	liquid (Supplied in 10 mM Tris-HCl pH 7.4, 400 mM KCl, 0.1 mM EDTA, 1 mM DTT, 200 µg/ml rAlbumin and 50% [v/v] glycerol)
Concentration	10 units/µl
Storage	store at -20 °C
Note	For research use only
Source	Diplococcus pneumoniae G41, recombinant, E.coli
Expiration Date	12 months from date of receipt.