

Product Datasheet

NaeI (orb93997)

Description	Isoschizomers: PdiI, Neoschizomers: MroNI, NgoMIV
Form/Appearance	liquid (Supplied in 10 mM Tris-HCl pH 7.4, 50 mM KCl, 0.1 mM EDTA, 1 mM DTT, 200 µg/ml rAlbumin and 50% [v/v] glycerol)
Concentration	10 units/µl
Storage	store at -20 °C
Note	For research use only
Source	Nocardia aerocolonigenes, recombinant, Streptomyces sp.
Expiration Date	12 months from date of receipt.