

Product Datasheet

BshFI (orb94016)

Description	Isoschizomers: BspANI, BsuRI, HaeIII, Pali, Phoi
Form/Appearance	liquid (Supplied in 10 mM Tris-HCl pH 7.4, 50 mM KCl, 0.1 mM EDTA, 1 mM DTT, 200 µg/ml rAlbumin and 50% [v/v] glycerol)
Concentration	10 units/µl
Storage	store at -20 °C
Note	For research use only
Source	Bacillus sphaericus
Expiration Date	12 months from date of receipt.