

Product Datasheet

Anti-FAM38A/PE-Cy5.5 (orb944465)

Description	Anti-FAM38A/PE-Cy5.5
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	PE/Cy5.5
Tested Applications	FC, ICC
Immunogen	KLH conjugated synthetic peptide derived from human FAM38A
Target	Piezo1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	287kDa
Uniprot ID	Q92508
Dilution Range	Flow-Cyt=1ug/Test,ICC=1:100
Expiration Date	12 months from date of receipt.