

Product Datasheet

Anti-GPCR ORF4/GPCR GPR137/Cy5 (orb944875)

Description	Anti-GPCR ORF4/GPCR GPR137/Cy5
Species/Host	Rabbit
Reactivity	Human
Conjugation	Cy5
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human GPCR ORF4/GPCR GPR137
Target	GPCR ORF4/GPCR GPR137
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	46kDa
Uniprot ID	Q96N19
Dilution Range	IF=1:100-500,ICC=1:100-500
Expiration Date	12 months from date of receipt.