


Product Datasheet

alpha 2a Adrenergic Receptor antibody [Out of stock] (orb96033)


Description	Rabbit polyclonal antibody against alpha 2a Adrenergic Receptor
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Immunogen	Synthetic peptide from human Adrenergic Receptor alpha-2A. (AA Range: 331-380)
Target	ADRA2A
Preservatives	PBS (without Mg ²⁺ , Ca ²⁺), pH 7.4, 150 mM NaCl, 0.02% Sodium Azide, 50% Glycerol
Concentration	1 mg/ml
Storage	Store at -20°C for up to 1 year.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Entrez	150
Dilution Range	ELISA:1:20000,IF/ICC:1:100-1:500,IHC-P:10 µg/ml,WB:1:500-1:1000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...



Immunofluorescence analysis of NIH-3T3 C...



Western blot analysis of Jurkat/LOVO cel...