

Product Datasheet

Anti-phospho-NMDAR2A (Ser1232)/PE-Cy5.5 [Out of stock] (orb960478)

Description	Anti-phospho-NMDAR2A (Ser1232)/PE-Cy5.5
Species/Host	Rabbit
Conjugation	PE/Cy5.5
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human NMDAR2A around the phosphorylation site of Ser1232
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	163kDa
Uniprot ID	Q00959 , P35436 , Q12879
Expiration Date	12 months from date of receipt.