


Product Datasheet

Luteinizing Hormone antibody [Out of stock] (orb96101)

Description	Mouse Monoclonal to Luteinizing Hormone.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Purified human LH-Beta.
Target	LHB
Preservatives	10 mM PBS, pH 7.4, 0.09% Sodium Azide, 0.2% rAlbumin
Concentration	0.2 mg/ml
Storage	Store at 2 to 8°C. Stable for 24 months.
Note	For research use only
Isotype	IgG1,k
Clonality	Monoclonal
Entrez	3972
Dilution Range	IHC-P:20 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...