

Product Datasheet

Anti-phospho-STAT3 (Ser727)/PE-Cy5 [Out of stock] (orb961140)

Description	Anti-phospho-STAT3 (Ser727)/PE-Cy5 [Out of stock]
Species/Host	Rabbit
Conjugation	Cy5
Tested Applications	FC, ICC, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human STAT3 around the phosphorylation site of Ser727
Target	phospho-STAT3 (Ser727)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P40763
Dilution Range	WB: 1:500-2000, IHC-P: 1:100-200, IHC-F: 1:100-200, ICC: 1:50-100, Flow-Cyt: 2ug/Test
Expiration Date	12 months from date of receipt.