

Product Datasheet

Anti-phospho-Smad2 (Ser467)/Cy5 (orb962071)

Description	Anti-phospho-Smad2 (Ser467)/Cy5
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Cy5
Tested Applications	ICC
Immunogen	KLH conjugated synthesised phosphopeptide derived from human Smad2 around the phosphorylation site of Ser467
Target	phospho-Smad2 (Ser467)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	58kDa
Uniprot ID	Q15796
Expiration Date	12 months from date of receipt.