

Product Datasheet

Anti-phospho-PDCD4 (Ser67)/AP (orb962895)

Product Description

Anti-phospho-PDCD4 (Ser67)/AP

Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	AP
Tested Applications	ICC
Immunogen	KLH conjugated synthesised phosphopeptide derived from human PDCD4 around the phosphorylation site of Ser67
Target	PDCD4
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	51kDa
Uniprot ID	Q53EL6
Dilution Range	WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Expiration Date	12 months from date of receipt.