

Product Datasheet

Anti-Phospho-Insulin Receptor Beta (Tyr1361)/Cy5 (orb964698)

Description	Anti-Phospho-Insulin Receptor Beta (Tyr1361)/Cy5
Species/Host	Rabbit
Reactivity	Human
Conjugation	Cy5
Tested Applications	FC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Insulin Receptor Beta around the phosphorylation site of Tyr1361
Target	Phospho-Insulin Receptor Beta (Tyr1361)
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
CAS Number	3
Uniprot ID	P06213
Expiration Date	12 months from date of receipt.