

Product Datasheet

Anti-Phospho-Insulin Receptor Beta (Tyr1361)/PE-Cy5.5 (orb965907)

Description

Anti-Phospho-Insulin Receptor Beta

Species/Host	Rabbit
Reactivity	Human
Conjugation	PE/Cy5.5
Tested Applications	FC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Insulin Receptor Beta around the phosphorylation site of Tyr1361
Target	Phospho-Insulin Receptor Beta (Tyr1361)
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
CAS Number	3
Uniprot ID	P06213
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558