

Product Datasheet

Anti-Anti-phospho-Smad2 (Ser250)/AP (orb966210)

Description	Anti-Anti-phospho-Smad2 (Ser250)/AP
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	AP
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human Smad2 around the phosphorylation site of Ser250
Target	phospho-Smad2 (Ser250)
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q15796
Dilution Range	WB=1:500-2000,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.