

Product Datasheet

Anti-Smac/AP (orb966317)

Description	Anti-Smac/AP
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	AP
Tested Applications	ICC, WB
Immunogen	KLH conjugated synthetic peptide derived from human Smac
Target	Smac
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	21kDa
Uniprot ID	Q9NR28
Dilution Range	WB=1:500-2000,ICC=1:100
Expiration Date	12 months from date of receipt.