

Product Datasheet

Anti-Anti-phospho-Smad2 (Ser250)/PE-Cy7 (orb967479)

Description

Anti-Anti-phospho-Smad2

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	PE/Cy7
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human Anti-phospho-Smad2 (Ser250)
Target	SMAD2
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	58kDa
Uniprot ID	Q15796
Expiration Date	12 months from date of receipt.