

Product Datasheet

Anti-caspase-3 p17 subunit /PE-Cy7 (orb967698)

Description	Anti-caspase-3 p17 subunit /PE-Cy7
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	PE/Cy7
Tested Applications	FC, ICC
Immunogen	KLH conjugated synthetic peptide derived from human caspase-3 p17 subunit
Target	caspase-3 p17 subunit
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P42574
Dilution Range	Flow-Cyt=1µg/Test, ICC=1:100-500
Expiration Date	12 months from date of receipt.