


Product Datasheet

Human CXCL6/GCP2 ELISA Kit (orb96772)

Description	Human CXCL6/GCP2 ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: V40-N114
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	11897 MW
Uniprot ID	P80162
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human
CXCL6/GCP2
PicoKine
ELISA Kit
stan...