


Product Datasheet

Rat IFN Gamma/IFNG/Interferon Gamma ELISA Kit (orb96773)

Description	Rat IFN Gamma/IFNG/Interferon Gamma ELISA Kit
Reactivity	Rat
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: E23-C156
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants and serum.
MW	17918 MW
Uniprot ID	P01581
Sensitivity	6 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum is less than the lowest standard, 31.2 pg/ml.
Expiration Date	Please enquire.



Rat IFN
gamma
PicoKine
ELISA Kit
standar...