


Product Datasheet

Human OSM/Oncostatin M ELISA Kit (orb96776)

Description	Human OSM/Oncostatin M ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: A26R-220
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	28484 MW
Uniprot ID	P13725
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 20-300 pg/ml.
Expiration Date	Please enquire.



Human
OSM/Oncostatin
M PicoKine
ELISA Ki...