


Product Datasheet

Mouse ACE/Cd143 ELISA Kit (orb96777)

Description	Mouse ACE/Cd143 ELISA Kit
Reactivity	Mouse
Range	93.8 pg/ml - 6,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: L35-Q1264
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin).
MW	150918 MW
Uniprot ID	P09470
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:500, concentration in serum and plasma is around 500 ng/ml.
Expiration Date	Please enquire.



Mouse ACE
PicoKine
ELISA Kit
standard
cu...