

Product Datasheet

Anti-Synaptopodin/PE-Cy5.5 (orb967796)

Description	Anti-Synaptopodin/PE-Cy5.5
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	PE/Cy5.5
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human Synaptopodin
Target	Synaptopodin
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8N3V7
Expiration Date	12 months from date of receipt.