


Product Datasheet

Mouse CTLA4 ELISA Kit (orb96781)

Description	Mouse CTLA4 ELISA Kit
Reactivity	Mouse
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A37-F162
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin).
MW	24993 MW
Uniprot ID	P09793
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 0-200 pg/ml.
Expiration Date	Please enquire.



Mouse
CTLA4
PicoKine
ELISA Kit
standard ...