


Product Datasheet

Human CXCL5/LIX/ENA-78 ELISA Kit (orb96782)

Description	Human CXCL5/LIX/ENA-78 ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: A37-N114
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
MW	11972 MW
Uniprot ID	P42830
Sensitivity	10 pg/ml
Dilution Range	For serum, dilution ratio is 1:1 and concentration is 442-3544 pg/ml. For plasma, dilution ratio is 1:1 and concentration is 0-1082 pg/ml.
Expiration Date	Please enquire.



Human CXCL5/ENA-78 PicoKine ELISA Kit st...