


Product Datasheet

Human CXCL16 ELISA Kit (orb96783)

Description	Human CXCL16 ELISA Kit
Reactivity	Human
Range	93.7 pg/ml - 6,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: N49-P137
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	27579 MW
Uniprot ID	Q9H2A7
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:5, concentration in serum and plasma is around 2 ng/ml.
Expiration Date	Please enquire.



Human
CXCL16
PicoKine
ELISA Kit
standard...