


Product Datasheet

Mouse Decorin ELISA Kit (orb96784)

Description	Mouse Decorin ELISA Kit
Reactivity	Mouse
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: D31-K354
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	39809 MW
Uniprot ID	P28654
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Mouse
Decorin
PicoKine
ELISA Kit
standar...