


Product Datasheet

Human HBEGF ELISA Kit (orb96788)

Description	Human HBEGF ELISA Kit
Reactivity	Human
Range	7.8 pg/ml - 500 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: D63-L148
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	culture supernatants, cell lysates, serum, plasma (heparin, EDTA), saliva, urine and human milk.
MW	23067 MW
Uniprot ID	Q99075
Sensitivity	10 pg/ml
Dilution Range	For serum, dilution ratio is 1:1 and concentration is 36-451 pg/ml. For plasma, dilution ratio is 1:1 and concentration is 0-114 pg/ml.
Expiration Date	Please enquire.



Human
HBEGF
PicoKine
ELISA Kit
standard ...