


Product Datasheet

Human IL-6R Alpha / IL 6 Receptor ELISA Kit (orb96789)

Description	Human IL-6R Alpha / IL 6 Receptor ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: L20-D358
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA, citrate) and urine.
MW	51548 MW
Uniprot ID	A0N0L5
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1000, concentration in serum and plasma is 10837-53930 pg/ml.
Expiration Date	Please enquire.



Human IL-6R alpha PicoKine ELISA Kit sta...