


Product Datasheet

Human CSF1R/M-CSFR ELISA Kit (orb96790)

Description	Human CSF1R/M-CSFR ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: I20-E512
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	107984 MW
Uniprot ID	P07333
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:50, concentration in serum and plasma is around 65 ng/ml.
Expiration Date	Please enquire.



Human
CSF1R/M-
CSFR
PicoKine
ELISA Kit
st...