


Product Datasheet

Mouse CSF1R/M-CSFR ELISA Kit (orb96791)

Description	Mouse CSF1R/M-CSFR ELISA Kit
Reactivity	Mouse
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A20-S511
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	109179 MW
Uniprot ID	P09581
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1000, concentration in serum and plasma is 400-947 ng/ml.
Expiration Date	Please enquire.



Mouse
CSF1R/M-
CSFR
PicoKine
ELISA Kit
st...