


Product Datasheet

Human LOX-1/OLR1 ELISA Kit (orb96792)

Description	Human LOX-1/OLR1 ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: S61-Q273
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	30959 MW
Uniprot ID	P78380
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 300 pg/ml.
Expiration Date	Please enquire.



Human LOX-1/OLR1 PicoKine ELISA Kit stan...