


Product Datasheet

Mouse CD10/Neprilysin ELISA Kit (orb96795)

Description	Mouse CD10/Neprilysin ELISA Kit
Reactivity	Mouse
Range	93.7 pg/ml - 6,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Y52-W750
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	85702 MW
Uniprot ID	Q61391
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 200 pg/ml.
Expiration Date	Please enquire.



Mouse
Neprilysin/CD10
PicoKine ELISA
Kit...