


Product Datasheet

Human PTX3/Pentraxin 3 ELISA Kit (orb96796)

Description	Human PTX3/Pentraxin 3 ELISA Kit
Reactivity	Human
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: E18-S381
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (EDTA) and saliva.
MW	41976 MW
Uniprot ID	P26022
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 350 pg/ml.
Expiration Date	Please enquire.



Human
PTX3/Pentraxin
3 PicoKine
ELISA Ki...