


Product Datasheet

Mouse PTX3/Pentraxin 3 ELISA Kit (orb96797)

Description	Mouse PTX3/Pentraxin 3 ELISA Kit
Reactivity	Mouse
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: E18-S381
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (EDTA).
MW	41812 MW
Uniprot ID	P48759
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:20, concentration in serum and plasma is around 36 ng/ml.
Expiration Date	Please enquire.



Mouse
PTX3/Pentraxin
3 PicoKine
ELISA Ki...