

Product Datasheet

Anti-LC3A/B/PE-Cy5.5 (orb968014)

Description	Anti-LC3A/B/PE-Cy5.5
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	PE/Cy5.5
Tested Applications	FC, ICC
Immunogen	KLH conjugated synthetic peptide derived from human LC3A/B
Target	LC3A/B
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9GZQ8
Dilution Range	Flow-Cyt=1µg/Test, ICC=1:100
Expiration Date	12 months from date of receipt.