


Product Datasheet

Human CD244/2B4 ELISA Kit (orb96808)

Description	Human CD244/2B4 ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: C22-R221
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	41616 MW
Uniprot ID	Q9BZW8
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 15.6 pg/ml.
Expiration Date	Please enquire.



Human
CD244/2B4
PicoKine
ELISA Kit
stand...