

Product Datasheet

HIF-2-alpha/EPAS1 Antibody (orb96817)

Description

HIF-2-alpha/EPAS1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

A synthetic peptide corresponding to a sequence in the middle region of human HIF-2-alpha (282-301aa, HALDSENMTKSHQNLCTKGQ), identical to the related mouse and rat sequence.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Predicted Species: Species predicted to be fit for the product based on sequence similarities. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW


96718 MW

Uniprot ID
Q9JHS1
Dilution Range


Western blot, 0.1-0.5µg/ml, Human, Mouse, Rat

Expiration Date


12 months from date of receipt.



Western blot analysis of HIF2A using ant...



Western blot analysis of HIF-2A using an...



Western blot analysis of HIF-2A using an...