

Product Datasheet

Anti-Cleaved PARP/AP (orb968439)

Description	Anti-Cleaved PARP/AP
Species/Host	Mouse
Reactivity	Human
Conjugation	AP
Tested Applications	ELISA, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human PARP
Target	Cleaved PARP1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P09874
Dilution Range	IHC-P=1:100-500,IHC-F=1:500-1000,ELISA=1:200-10000
Expiration Date	12 months from date of receipt.