

Product Datasheet

Anti-HP1 gamma/AP (orb968810)

Description	Anti-HP1 gamma/AP
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	AP
Tested Applications	ELISA, ICC, IHC-Fr, IHC-P, WB
Immunogen	Recombinant human HP1 gamma Protein
Target	HP1 gamma
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	Q13185
Dilution Range	WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,ICC=1:100-500,ELISA=1:1000-5000
Expiration Date	12 months from date of receipt.