

Product Datasheet

Anti-Ki-67/PE-Cy5.5 (orb969320)

Description	Anti-Ki-67/PE-Cy5.5
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	PE/Cy5.5
Tested Applications	FC, ICC
Immunogen	KLH conjugated synthetic peptide derived from mouse Ki-67
Target	Ki-67
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P46013
Dilution Range	Flow-Cyt=1ug/Test,ICC=1:100
Expiration Date	12 months from date of receipt.