

Product Datasheet

Anti-phospho-p38 MAPK (Tyr323)/Cy5 [Out of stock] (orb969524)

Description

Anti-phospho-p38 MAPK (Tyr323)/Cy5

Species/Host

Rabbit

Conjugation

Cy5

Tested Applications

ICC, IF

Immunogen

KLH conjugated synthetic peptide derived from human phospho-p38 MAPK (Tyr323)

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

1mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

41kDa

Uniprot ID
Q16539, P70618, P47811, O02812
Dilution Range

ICC=1:50-200 IF=1:50-200 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.

Expiration Date

12 months from date of receipt.