

Product Datasheet

BCRP/ABCG2 Antibody (orb96996)

Description

BCRP/ABCG2 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

A synthetic peptide corresponding to a sequence in the middle region of human BCRP/ABCG2 (154-168aa NHEKNERINRVIQEL).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

WB: The detection limit for ABCG2 is approximately 1ng/lane under reducing conditions. Tested Species: In-house tested species with positive results. Predicted Species: Species predicted to be fit for the product based on sequence similarities. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

72314 MW

Uniprot ID

[Q9UNQ0](#)


Dilution Range

Western blot, 0.1-0.5µg/ml, Human

Expiration Date

12 months from date of receipt.

1 2 3 4 5
— 250KD
— 130KD



Western blot analysis of
Lane 1: HELA
Ce...