


Product Datasheet

BCRP/ABCG2 Antibody (orb96997)

Description

BCRP/ABCG2 Antibody

Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide corresponding to a sequence in the middle region of mouse BCRP/ABCG2 (152-167aa KNHEKNERINTIIKEL), different from the related rat sequence by one amino acid.
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	WB: The detection limit for ABCG2 is approximately 1ng/lane under reducing conditions. Tested Species: In-house tested species with positive results. Predicted Species: Species predicted to be fit for the product based on sequence similarities. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	72978 MW
Uniprot ID	Q7TMS5
Dilution Range	Western blot, 0.1-0.5µg/ml, Mouse, Rat
Expiration Date	12 months from date of receipt.



Western blot analysis of
Lane 1: Rat
Liv...