

Product Datasheet

Anti-CD41/Integrin alpha 2b/PE-Cy5.5 (orb970260)

Description	Anti-CD41/Integrin alpha 2b/PE-Cy5.5
Species/Host	Mouse
Reactivity	Human
Conjugation	PE/Cy5.5
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human ITGA2B
Target	CD41/ITGA2B
Preservatives	0.01M PBS(pH7.4) with 0.1% Proclin300
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P08514
Expiration Date	12 months from date of receipt.