

Product Datasheet

MEKK3/MAP3K3 Antibody (orb97050)

Description

MEKK3/MAP3K3 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

A synthetic peptide corresponding to a sequence at the N-terminus of human MEKK3 (11-25aa MNDLVALQMNRHRM).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

WB: The detection limit for MAP3K3 is approximately 0.5ng/lane under reducing conditions. Tested Species: In-house tested species with positive results. Predicted Species: Species predicted to be fit for the product based on sequence similarities. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

70898 MW

Uniprot ID

[Q99759](https://www.uniprot.org/entry/Q99759)

Dilution Range

Western blot, 0.1-0.5µg/ml, Human

Expiration Date

12 months from date of receipt.

1 2 3 4 — 250KD
— 130KD

— 70KD
— 55KD
— 35KD
— 25KD
— 15KD

Western blot
analysis of
Lane 1: HELA
Ce...