

Product Datasheet

Anti-phospho-p38 MAPK (Tyr323)/PE-Cy5.5 [Out of stock] (orb970737)

Description

Anti-phospho-p38 MAPK (Tyr323)/PE-

Species/Host	Rabbit
Conjugation	PE/Cy5.5
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human phospho-p38 MAPK (Tyr323)
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.05M PB, pH 7.5.
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	41kDa
Uniprot ID	Q16539 , P70618 , P47811 , O02812
Dilution Range	ICC=1:50-200 IF=1:50-200 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Expiration Date	12 months from date of receipt.