

Product Datasheet

Anti-phospho-Gli1 (S112/ T115/ S116)/AP (orb971039)

Description	Anti-phospho-Gli1 (S112/ T115/ S116)/AP
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	AP
Tested Applications	ELISA, ICC, IHC-Fr, IHC-P
Immunogen	KLH conjugated sSynthesised phosphopeptidepeptide derived from human Gli1 around the phosphorylation site of S112/ T115/ S116
Target	phospho-Gli1 (S112/ T115/ S116)
Preservatives	PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1mg/1ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P08151
Dilution Range	IHC-P=1:100-500,IHC-F=1:100-500,ICC=1:100-500,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.