

Product Datasheet

Anti-SCARB1/Scavenger Receptor BI /PE-Cy5.5 [Out of stock] (orb971682)

Description	Anti-SCARB1/Scavenger Receptor BI /PE-Cy5.5
Species/Host	Rabbit
Conjugation	PE/Cy5.5
Immunogen	KLH conjugated synthetic peptide derived from mouse SCARB1/Scavenger Receptor BI
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.05M PB, pH 7.5.
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	61kDa
Uniprot ID	P97943, Q61009, Q8WTV0
Expiration Date	12 months from date of receipt.